

## **Product datasheet for TA505894**

### **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

## **HGS Mouse Monoclonal Antibody [Clone ID: OTI3H7]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3H7
Applications: IF, WB

**Reactivity:** WB 1:4000, IF 1:100 **Reactivity:** Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HGS(NP\_004703) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 86 kDa

**Gene Name:** hepatocyte growth factor-regulated tyrosine kinase substrate

Database Link: NP 004703

Entrez Gene 15239 MouseEntrez Gene 56084 RatEntrez Gene 9146 Human

<u>014964</u>

**Background:** The protein encoded by this gene regulates endosomal sorting and plays a critical role in the

recycling and degradation of membrane receptors. The encoded protein sorts

monoubiquitinated membrane proteins into the multivesicular body, targeting these proteins

for lysosome-dependent degradation. [provided by RefSeq, Dec 2010]

Synonyms: HRS







**Protein Families:** Druggable Genome

**Protein Pathways:** Endocytosis

# **Product images:**





HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HGS (Cat# [RC200609], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HGS(Cat# TA505894). Positive lysates [LY401488] (100ug) and [LC401488] (20ug) can be purchased separately from OriGene.

Anti-HGS mouse monoclonal antibody (TA505894) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HGS ([RC200609]).