

Product datasheet for TA506201S

ERG Mouse Monoclonal Antibody [Clone ID: OTI8A9]

Product data:

Product Type: Primary Antibodies Clone Name: OTI8A9 IF, IHC, WB **Applications:** Recommended Dilution: WB: 1:4000, IHC: 1:150, IF: 1:50-1:100 **Reactivity:** Human, Mouse, Rat Host: Mouse Isotype: lgG2b **Clonality:** Monoclonal Full length human recombinant protein of human ERG(NP_891548) produced in HEK293T Immunogen: cell. Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. Concentration: 1 mg/ml **Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) **Conjugation:** Unconjugated Store at -20°C as received. Storage: Stability: Stable for 12 months from date of receipt. Predicted Protein Size: 53.7 kDa Gene Name: ETS transcription factor ERG Database Link: NP 891548 Entrez Gene 13876 MouseEntrez Gene 170909 RatEntrez Gene 2078 Human P11308



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Serigene ERG Mouse Monoclonal Antibody [Clone ID: OTI8A9] – TA506201S

Background:This gene encodes a member of the erythroblast transformation-specific (ETS) family of
transcriptions factors. All members of this family are key regulators of embryonic
development, cell proliferation, differentiation, angiogenesis, inflammation, and apoptosis.
The protein encoded by this gene is mainly expressed in the nucleus. It contains an ETS DNA-
binding domain and a PNT (pointed) domain which is implicated in the self-association of
chimeric oncoproteins. This protein is required for platelet adhesion to the subendothelium,
inducing vascular cell remodeling. It also regulates hematopoesis, and the differentiation and
maturation of megakaryocytic cells. This gene is involved in chromosomal translocations,
resulting in different fusion gene products, such as TMPSSR2-ERG and NDRG1-ERG in
prostate cancer, EWS-ERG in Ewing's sarcoma and FUS-ERG in acute myeloid leukemia.
Multiple alternatively spliced transcript variants encoding different isoforms have been
identified. [provided by RefSeq, Jan


Synonyms:

erg-3; p55


Protein Families:

Druggable Genome, Transcription Factors

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ERG ([RC208093], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ERG. Positive lysates [LY403650] (100ug) and [LC403650] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-ERG monoclonal antibody (1:200).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US








Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-ERG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA506201])

Anti-ERG mouse monoclonal antibody ([TA506201]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ERG ([RC208093]) at 1:100

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Immunofluorescent staining of HeLa cells using anti-ERG mouse monoclonal antibody ([TA506201]) at 1:100

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US