

# **Product datasheet for TA506343AM**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CD64 (FCGR1A) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI3D3]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3D3

**Applications:** FC, IF, IHC, WB

Recommended Dilution: WB 1:4000, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human FCGR1A(NP\_000557) produced in HEK293T

cell.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Conjugation:** Biotin

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 42.5 kDa

**Gene Name:** Fc fragment of lgG receptor la

Database Link: NP 000557

Entrez Gene 2209 Human

P12314

**Background:** This gene encodes a protein that plays an important role in the immune response. This

protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family

members located on chromosome 1. [provided by RefSeq, Jul 2008]

Synonyms: CD64; CD64A; FCRI; IGFR1

**Protein Families:** Transmembrane


Protein Pathways: Fc gamma R-mediated phagocytosis, Hematopoietic cell lineage, Systemic lupus

erythematosus







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY FCGR1A (Cat# [RC207487], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FCGR1A(Cat# [TA506343]). Positive lysates [LY400193] (100ug) and [LC400193] (20ug) can be purchased separately from OriGene.



Anti-FCGR1A mouse monoclonal antibody ([TA506343]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FCGR1A ([RC207487]).



Flow cytometric analysis of living 293T cells transfected with FCGR1A overexpression plasmid ([RC207487]), Red)/empty vector ([PS100001], Blue) using anti-FCGR1A antibody ([TA506343]). Cells incubated with a non-specific antibody (Green) were used as isotype control (1:100).