

Product datasheet for TA506590

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

beta II Tubulin (TUBB2A) Mouse Monoclonal Antibody [Clone ID: OTI3E6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3E6

Applications: IF, IHC, WB

Recommended Dilution: WB 1:4000, IF 1:100, IHC: 1:150

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TUBB2A(NP_001060) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: tubulin beta 2A class IIa

Database Link: NP 001060

Entrez Gene 22151 MouseEntrez Gene 498736 RatEntrez Gene 7280 Human

Q13885

Synonyms: CDCBM5; dJ40E16.7; TUBB; TUBB2


Protein Families: Druggable Genome

Protein Pathways: Gap junction, Pathogenic Escherichia coli infection







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TUBB2A (Cat# [RC205881], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TUBB2A(Cat# TA506590). Positive lysates [LY420706] (100ug) and [LC420706] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-TUBB2A mouse monoclonal antibody. (TA506590)



Anti-TUBB2A mouse monoclonal antibody (TA506590) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TUBB2A ([RC205881]).





Immunofluorescent staining of HeLa cells using anti-TUBB2A mouse monoclonal antibody (TA506590).