

Product datasheet for **TA507002AM**


MFAP3 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1D3]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1D3
Applications:	IF, WB
Recommended Dilution:	WB 1:4000, IF 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human MFAP3(NP_005918) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	40 kDa
Gene Name:	microfibril associated protein 3
Database Link:	NP_005918 Entrez Gene 4238 Human P55082
Synonyms:	DKFZp586F0219
Protein Families:	Transmembrane



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MFAP3 ([RC204656], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MFAP3. Positive lysates [LY416962] (100ug) and [LC416962] (20ug) can be purchased separately from OriGene.



Anti-MFAP3 mouse monoclonal antibody ([TA507002]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MFAP3 ([RC204656]).