

Product datasheet for TA507028BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Hydroxysteroid (17 beta) Dehydrogenase 4 (HSD17B4) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4C4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4C4

Applications: IF, IHC, WB

Recommended Dilution: WB 1:4000, IF 1:100, IHC 1:150

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HSD17B4(NP_000405) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: hydroxysteroid 17-beta dehydrogenase 4

Database Link: NP 000405

Entrez Gene 15488 MouseEntrez Gene 3295 Human

P51659

Background: The protein encoded by this gene is a bifunctional enzyme that is involved in the peroxisomal

beta-oxidation pathway for fatty acids. It also acts as a catalyst for the formation of 3-ketoacyl-CoA intermediates from both straight-chain and 2-methyl-branched-chain fatty acids. Defects in this gene that affect the peroxisomal fatty acid beta-oxidation activity are a cause of D-bifunctional protein deficiency (DBPD). An apparent pseudogene of this gene is

present on chromosome 8. [provided by RefSeq, Jul 2008]






Hydroxysteroid (17 beta) Dehydrogenase 4 (HSD17B4) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4C4] – TA507028BM

Synonyms: DBP; MFE-2; MPF-2; PRLTS1; SDR8C1


Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Primary bile acid biosynthesis

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HSD17B4 (Cat# [RC200460], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HSD17B4(Cat# [TA507028]). Positive lysates [LY424737] (100ug) and [LC424737] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.





Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human Ovary tissue within the normal limits using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.






Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:150 dilution.





Anti-HSD17B4 mouse monoclonal antibody ([TA507028]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HSD17B4 ([RC200460]).



Immunofluorescent staining of HeLa cells using anti-HSD17B4 mouse monoclonal antibody ([TA507028]) at 1:100 dilution.