

Product datasheet for **TA507060AM**

UCHL5IP (HAUS7) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1E8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1E8
Applications:	IF, IHC, WB
Recommended Dilution:	WB 1:4000, IHC 1:150, IF: 1:50-1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human HAUS7(NP_059988) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	HAUS augmin like complex subunit 7
Database Link:	NP_059988 Entrez Gene 55559 Human Q99871
Synonyms:	UCHL5IP; UIP1




[View online »](#)


Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HAUS7 (Cat# [RC203737], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HAUS7(Cat# [TA507060]). Positive lysates [LY413731] (100ug) and [LC413731] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffin-embedded Human breast tissue within the normal limits using anti-HAUS7 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA507060])



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-HAUS7 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA507060])



Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-HAUS7 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA507060])



Immunofluorescent staining of HeLa cells using anti-HAUS7 mouse monoclonal antibody ([TA507060]) at 1:100