

Product datasheet for **TA507110BM**

Angiomotin (AMOT) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1F11]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1F11
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human AMOT(NP_573572) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	72.4 kDa
Gene Name:	angiomotin
Database Link:	NP_573572 Entrez Gene 27494 Mouse Entrez Gene 300289 Rat Entrez Gene 154796 Human Q4VCS5
Synonyms:	KIAA1071



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AMOT ([RC220347], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AMOT. Positive lysates [LY408854] (100ug) and [LC408854] (20ug) can be purchased separately from OriGene.