

## Product datasheet for **TA507207S**

### **PAAF1 Mouse Monoclonal Antibody [Clone ID: OTI1F11]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI1F11
Applications:	IF, WB
Recommended Dilution:	WB 1:400~4000, IF 1:100
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human PAAF1(NP_079431) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42 kDa
Gene Name:	proteasomal ATPase associated factor 1
Database Link:	<a href="#">NP_079431</a> <a href="#">Entrez Gene 80227 Human</a> <a href="#">Q9BRP4</a>
Synonyms:	PAAF; Rpn14; WDR71




[View online »](#)

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PAAF1 (Cat# [RC202994], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PAAF1(Cat# [TA507207]). Positive lysates [LY410872] (100ug) and [LC410872] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-PAAF1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Anti-PAAF1 mouse monoclonal antibody ([TA507207]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PAAF1 ([RC202994]).