

Product datasheet for **TA507234**

MINPP1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1B2
Applications:	IF, IHC, WB
Recommended Dilution:	WB 1:4000, IF 1:100, IHC 1:150
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human MINPP1(NP_004888) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54.9 kDa
Gene Name:	multiple inositol-polyphosphate phosphatase 1
Database Link:	NP_004888 Entrez Gene 9562 Human Q9UNW1
Background:	This gene encodes multiple inositol polyphosphate phosphatase; an enzyme that removes 3-phosphate from inositol phosphate substrates. It is the only enzyme known to hydrolyze inositol pentakisphosphate and inositol hexakisphosphate. This enzyme also converts 2,3 bisphosphoglycerate (2,3-BPG) to 2-phosphoglycerate; an activity formerly thought to be exclusive to 2,3-BPG synthase/2-phosphatase (BPGM) in the Rapoport-Luebering shunt of the glycolytic pathway. [provided by RefSeq, Sep 2009]




[View online »](#)


Synonyms: HIPER1; MINPP2; MIPP
Protein Families: Druggable Genome
Protein Pathways: Inositol phosphate metabolism

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY MINPP1 (Cat# [RC207581], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MINPP1 (Cat# TA507234). Positive lysates [LY417667] (100ug) and [LC417667] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-MINPP1 mouse monoclonal antibody (TA507234) at 1:150 dilution. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.



Anti-MINPP1 mouse monoclonal antibody (TA507234) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MINPP1 ([RC207581]).