

## Product datasheet for **TA507285AM**


### HSCB Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2B7]

#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Clone Name:             | OTI2B7  |
| Applications:           | WB  |
| Recommended Dilution:   | WB 1:4000   |
| Reactivity:             | Human   |
| Host:                   | Mouse   |
| Isotype:                | IgG2a   |
| Clonality:              | Monoclonal  |
| Immunogen:              | Full length human recombinant protein of human HSCB(NP_741999) produced in HEK293T cell                   |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                      |
| Concentration:          | 0.5 mg/ml   |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation:            | Biotin  |
| Storage:                | Store at -20°C as received.   |
| Stability:              | Stable for 12 months from date of receipt.  |
| Predicted Protein Size: | 27.2 kDa  |
| Gene Name:              | HscB mitochondrial iron-sulfur cluster cochaperone  |
| Database Link:          | <a href="#">NP_741999</a><br><a href="#">Entrez Gene 150274 Human</a><br><a href="#">Q8IWL3</a>           |
| Synonyms:               | DNAJC20; HSC20; JAC1  |



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HSCB ([RC200843], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HSCB. Positive lysates [LY406803] (100ug) and [LC406803] (20ug) can be purchased separately from OriGene.