

## Product datasheet for **TA507304AM**


### **ANKMY2 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2C3]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI2C3
Applications:	WB
Recommended Dilution:	WB 1:400~4000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ANKMY2(NP_064715) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	49.1 kDa
Gene Name:	ankyrin repeat and MYND domain containing 2
Database Link:	<a href="#">NP_064715</a> <a href="#">Entrez Gene 217473 Mouse</a> <a href="#">Entrez Gene 314046 Rat</a> <a href="#">Entrez Gene 57037 Human</a> <a href="#">Q8IV38</a>
Synonyms:	ZMYND20



[View online »](#)

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ANKMY2 ([RC206770], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ANKMY2. Positive lysates [LY412559] (100ug) and [LC412559] (20ug) can be purchased separately from OriGene.