

Product datasheet for TA590105

NSE (ENO2) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, WB

Recommended Dilution: WB 1:20000, ELISA 1:100-1:2000

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: DNA immunization. This antibody is specific for the N Terminus Region of the target protein.

Formulation: 20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0

Concentration: 0.90292 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: enolase 2

Database Link: NP 001966

Entrez Gene 2026 Human

P09104

Background: Neuron-specific enolase (NSE) is a glycolytic isoenzyme which is located in central and

peripheral neurons and neuroendocrine cells. This enzyme is released into the CSF when neural tissue is injured. Neoplasms derived from neural or neuroendocrine tissue may release NSE into the blood. NSE is a useful substance that has been detected in patients with certain tumors, namely: neuroblastoma, small cell lung cancer, medullary thyroid cancer, carcinoid tumors, pancreatic endocrine tumors, and melanoma. Recombinant NSE was

expressed in E.coli and purified by conventional chromatography techniques.

Synonyms: HEL-S-279; NSE



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Note:


This antibody was generated by SDIX's Genomic Antibody Technology ® (GAT). Learn about

<u>GAT</u>


Protein Pathways:

Glycolysis / Gluconeogenesis, Metabolic pathways, RNA degradation

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ENO2 ([RC201085], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ENO2. Positive lysates [LY400724] (100ug) and [LC400724] (20ug) can be purchased separately from OriGene.



ELISA: NSE Antibody