

## Product datasheet for **TA591007S**


### Rabbit monoclonal anti-Dendra2 Antibody, clone OTIR5C7

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Rabbit monoclonal anti-Dendra2 Antibody, clone OTIR5C7
Clone Name:	OR-5C7
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Dendra2
Host:	Rabbit
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of Dendra2 expressed in HEK293 cells
Isotype:	IgG
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	Lot dependent; please refer to CoA along with shipment
Purification:	Purified from mouse ascites fluids or tissue/cell culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	27.5 kDa



[View online »](#)

**Product images:**

Western blot analysis of Dendra2 purified lysates (5ug) using anti-Dendra2 rabbit monoclonal antibody (TA591007) at 1:1000 dilution.