

## Product datasheet for **TA591008S**


### Rabbit monoclonal anti-His tag Antibody, clone OTIR5B12

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Rabbit monoclonal anti-His tag Antibody, clone OTIR5B12
Clone Name:	OR-5B12
Applications:	WB
Recommended Dilution:	WB: 1:500
Reactivity:	His tag
Host:	Rabbit
Clonality:	Monoclonal
Immunogen:	Synthetic peptide conjugated to KLH
Isotype:	IgG
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	Lot dependent; please refer to CoA along with shipment
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	ANTI HIS, ANTI HISTIDINE, ANTIHIS, anti-His, his tag antibody, Anti-His Antibody, his antibody, anti his tag, anti-his tag



[View online »](#)

**Product images:**

Western blot analysis of extracts from SF9 cell line (left lane) and SF9 cells transfected with C-terminal His-tagged gene (right lane), using anti-His tag rabbit monoclonal antibody ([TA591008]) at 1:500 dilution.