


Product datasheet for **TA591010**

Rabbit monoclonal anti-HA tag Antibody, clone OTIR5B11

Product data:

Product Type:	Tag Antibodies
Product Name:	Rabbit monoclonal anti-HA tag Antibody, clone OTIR5B11
Clone Name:	OTIR5B11
Applications:	WB
Recommended Dilution:	WB: 1:2000
Target:	Human
Host:	Rabbit
Clonality:	Monoclonal
Immunogen:	Synthetic peptide conjugated to KLH
Isotype:	IgG
Formulation:	1XPBS containing 50% glycerol, 1% BSA and 0.02% Sodium azide
Concentration:	0.5mg/mL
Conjugation:	Unconjugated
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.

Product images:



Western blot analysis of extracts from HEK293 cells transfected with C-terminal HA-tagged gene (lane 2, 3, 4, with 5, 10, 20uL from left to right) and negative control vector (lane 1) using anti-HA tag rabbit monoclonal antibody (TA591010) at 1:2000 dilution.



[View online »](#)