

## Product datasheet for **TA591029**

### **S100 alpha 2 (S100A2) Rabbit Monoclonal Antibody [Clone ID: OTIR4F6]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR4F6
Applications:	SISCAPA
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human S100A2 (NP_005969). The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	Lot dependent; please refer to CoA along with shipment
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	S100 calcium binding protein A2
Database Link:	<a href="#">NP_005969</a> <a href="#">Entrez Gene 6273 Human</a> <a href="#">P29034</a>
Synonyms:	CAN19; S100L



[View online »](#)