

Product datasheet for **TA591092**

CXCR4 Mouse Monoclonal Antibody [Clone ID: 772AA80]

Product data:


Product Type:	Primary Antibodies
Clone Name:	772AA80
Applications:	FC
Recommended Dilution:	Use at 0.1-1ug/ml. Flow
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	A near-full length cDNA gene (cytoplasmic C-terminus partially truncated) delivered by DNA immunization
Formulation:	PBS with 50% glycerol and 0.05% azide
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	C-X-C motif chemokine receptor 4
Database Link:	NP_003458 Entrez Gene 7852 Human P61073
Background:	SDF1 receptor
Synonyms:	CD184; D2S201E; FB22; HM89; HSY3RR; LAP-3; LAP3; LCR1; LESTR; NPY3R; NPYR; NPYRL; NPY3R; WHIM
Note:	Epitope N-terminus. This antibody was generated by SDIX's Genomic Antibody Technology® (GAT). Learn about GAT
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, GPCR, Transmembrane



[View online »](#)

Protein Pathways: Axon guidance, Chemokine signaling pathway, Cytokine-cytokine receptor interaction, Endocytosis, Leukocyte transendothelial migration


Product images:



Titration of CXCR4 antibody 772AA80 on Jurkat cells in flow cytometry

Clone name	Flow Cytometry			
	Jurkat	HEK 293	Ramos	CEM
AA80	++	+++	+++	+++

Relative staining of different cell types in flow cytometry



Epitope: N-terminus

Epitopes were defined by homologous scanning mutagenesis using human-mouse CXCR4 chimeras. Note: this method defines the key contact points and the complete set of contacts will likely be broader.