

Product datasheet for **TA592012**


P16(CDKN2A) Rabbit Monoclonal Antibody [Clone ID: OTIR5G10]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTIR5G10
Applications:	IHC
Recommended Dilution:	IHC 1:9000
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide corresponding to the C-terminus of Human P16(CDKN2A) (NP_000068) conjugated to KLH. The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.



[View online »](#)

Product images:

Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-P16 (CDKN2A) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 3 min, TA592012) (1:9000)



Immunohistochemical staining of paraffin-embedded Human cervical cancer tissue using anti-P16 (CDKN2A) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 3 min, TA592012) (1:9000)