

## Product datasheet for **TA592052**


### Desmin(DES) Rabbit Monoclonal Antibody [Clone ID: OTIR4A8]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTIR4A8
Applications:	IHC
Recommended Dilution:	IHC 1:2500
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide corresponding to the C-terminus of Human Desmin(DES) (NP_001918) conjugated to KLH. The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.




[View online »](#)

**Product images:**


Immunohistochemical staining of paraffin-embedded Human colon tissue within the normal limits using anti-Desmin (DES) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA592052) (1:2500)



Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-Desmin (DES) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA592052) (1:2500)



Immunohistochemical staining of paraffin-embedded Human gastric tissue within the normal limits using anti-Desmin (DES) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA592052) (1:2500)



Immunohistochemical staining of paraffin-embedded Human appendix tissue within the normal limits using anti-Desmin (DES) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, TA592052) (1:2500)