

## Product datasheet for TA592071S

### Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rockville, MD 20850, US

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

# AMH Rabbit Monoclonal Antibody [Clone ID: OTIR2H6]

## **Product data:**

**Product Type: Primary Antibodies** 

Clone Name: OTIR2H6

**Applications: ELISA** 

Recommended Dilution: ELISA 1:5000-10000

Reactivity: Human Host: Rabbit Isotype: lgG

Clonality: Monoclonal

Immunogen: Recombinant protein (26-560aa) of human AMH produced in HEK299 PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. Formulation:

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Predicted Protein Size:** 59.17 kDa

Gene Name: anti-Mullerian hormone

Database Link: NP 000470

Entrez Gene 268 Human

P03971



## AMH Rabbit Monoclonal Antibody [Clone ID: OTIR2H6] - TA592071S

**Background:** This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta)

superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene

expression. The encoded preproprotein is proteolytically processed to generate N- and C-terminal cleavage products that homodimerize and associate to form a biologically active noncovalent complex. This complex binds to the anti-Mullerian hormone receptor type 2 and

causes the regression of Mullerian ducts in the male embryo that would otherwise

differentiate into the uterus and fallopian tubes. This protein also plays a role in Leydig cell differentiation and function and follicular development in adult females. Mutations in this

gene result in persistent Mullerian duct syndrome. [provided by RefSeq, Jul 2016]

Synonyms: MIF; MIS

**Protein Families:** Druggable Genome, Secreted Protein

**Protein Pathways:** Cytokine-cytokine receptor interaction, TGF-beta signaling pathway