

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA592095S


# GS(GLUL) Rabbit Monoclonal Antibody [Clone ID: OTIR1A5]

## **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTIR1A5
Applications:	IHC
Recommended Dilution:	IHC 1:1500
Host:	Rabbit
lsotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide corresponding to the C-terminus of Human GS(GLUL) (NP_001028228) conjugated to KLH. The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.




### **Product images:**



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)




Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)

Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)


Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)



Immunohistochemical staining of paraffinembedded Human thyroid tissue within the normal limits using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)

Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592052]) (1:1500)

Immunohistochemical staining of paraffinembedded Human gastric tissue within the normal limits using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)



Immunohistochemical staining of paraffinembedded Human appendix tissue within the normal limits using anti-GS (GLUL) rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592095]) (1:1500)