

Product datasheet for TA592100S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

BNP (NPPB) Rabbit Monoclonal Antibody [Clone ID: OTIR1C4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTIR1C4

Applications: ELISA

Recommended Dilution: ELISA 1:5000-10000

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Monoclonal

Immunogen: Synthetic peptide (the amino acid sequence is considered to be commercially sensitive)

within Human BNP. The exact sequence is proprietary.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Predicted Protein Size: 14.7 kDa

Gene Name: natriuretic peptide B

Database Link: NP 002512

Entrez Gene 4879 Human

P16860

Synonyms: BNP; Iso-ANP

Protein Families: Druggable Genome, Secreted Protein, Stem cell - Pluripotency

