

## Product datasheet for **TA592106S**

### **BNP (NPPB) Rabbit Monoclonal Antibody [Clone ID: OTIR3G3]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR3G3
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human BNP. The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	14.7 kDa
Gene Name:	natriuretic peptide B
Database Link:	<a href="#">NP_002512</a> <a href="#">Entrez Gene 4879 Human</a> <a href="#">P16860</a>
Synonyms:	BNP; Iso-ANP
Protein Families:	Druggable Genome, Secreted Protein, Stem cell - Pluripotency



[View online »](#)