

Product datasheet for **TA592109**

BNP (NPPB) Rabbit Monoclonal Antibody [Clone ID: OTIR5G8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTIR5G8
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human BNP. The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	14.7 kDa
Gene Name:	natriuretic peptide B
Database Link:	NP_002512 Entrez Gene 4879 Human P16860
Synonyms:	BNP; Iso-ANP
Protein Families:	Druggable Genome, Secreted Protein, Stem cell - Pluripotency



[View online »](#)