

Product datasheet for **TA592120M**

MUC6 Rabbit Monoclonal Antibody [Clone ID: OTIR2F3]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTIR2F3
Applications:	IHC
Recommended Dilution:	IHC 1:4000
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human MUC6 (NP_005952). The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.




[View online »](#)

Product images:



Immunohistochemical staining of paraffin-embedded Human gastric tissue within the normal limits using anti-MUC6 rabbit monoclonal antibody. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.



Immunohistochemical staining of paraffin-embedded Human gastric carcinoma using anti-MUC6 rabbit monoclonal antibody. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.