

Product datasheet for **TA592122S**


MUC6 Rabbit Monoclonal Antibody [Clone ID: OTIR4C11]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTIR4C11
Applications:	IHC
Recommended Dilution:	IHC 1:500
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human MUC6 (NP_005952). The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.



[View online »](#)

Product images:

Immunohistochemical staining of paraffin-embedded Human gastric tissue within the normal limits using anti-MUC6 rabbit monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min, [TA592122]) (1:500)



Immunohistochemical staining of paraffin-embedded Human Gastric carcinoma using anti-MUC6 rabbit monoclonal antibody ([TA592122]). Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3 min. (1:500)