

## Product datasheet for **TA592127M**


### **HNF1-Beta(HNF1B) Rabbit Monoclonal Antibody [Clone ID: OTIR1E9]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR1E9
Applications:	IHC
Recommended Dilution:	IHC 1:500
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to the N-terminus of human HNF1-Beta(HNF1B) (NP_000449) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.



[View online »](#)

**Product images:**

Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-HNF1-Beta (HNF1B) rabbit monoclonal antibody. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.