

## Product datasheet for **TA592142**

### **gag Rabbit Monoclonal Antibody [Clone ID: OTIR1A4]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTIR1A4
<b>Applications:</b>	ELISA
<b>Recommended Dilution:</b>	ELISA 1:5000-10000
<b>Reactivity:</b>	HIV
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Recombinant protein fragment corresponding to amino acids 1-231 of HIV P24 produced in E.coli.
<b>Formulation:</b>	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	1 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Unconjugated
<b>Predicted Protein Size:</b>	25.58kDa



[View online »](#)