

## Product datasheet for **TA592146**

### **gag Rabbit Monoclonal Antibody [Clone ID: OTIR1E3]**

#### **Product data:**

|                                |   |
|--------------------------------|---|
| <b>Product Type:</b>           | Primary Antibodies  |
| <b>Clone Name:</b>             | OTIR1E3   |
| <b>Applications:</b>           | ELISA   |
| <b>Recommended Dilution:</b>   | ELISA 1:5000-10000  |
| <b>Reactivity:</b>             | HIV   |
| <b>Host:</b>                   | Rabbit  |
| <b>Isotype:</b>                | IgG   |
| <b>Clonality:</b>              | Monoclonal  |
| <b>Immunogen:</b>              | Recombinant protein fragment corresponding to amino acids 1-231 of HIV P24 produced in E.coli.            |
| <b>Formulation:</b>            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                      |
| <b>Concentration:</b>          | 1 mg/ml   |
| <b>Purification:</b>           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| <b>Conjugation:</b>            | Unconjugated  |
| <b>Predicted Protein Size:</b> | 25.58kDa  |



[View online »](#)