

Product datasheet for **TA592264**

gag Rabbit Monoclonal Antibody [Clone ID: OTIR2C5]

Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Clone Name: | OTIR2C5 |
| Applications: | ELISA |
| Recommended Dilution: | ELISA 1:5000-10000 |
| Reactivity: | HIV |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | Recombinant protein fragment corresponding to amino acids 1-231 of HIV P24 produced in E.coli. |
| Formulation: | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. |
| Concentration: | 1 mg/ml |
| Purification: | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation: | Unconjugated |
| Predicted Protein Size: | 25.58kDa |



[View online »](#)