

## Product datasheet for **TA592264S**

### **gag Rabbit Monoclonal Antibody [Clone ID: OTIR2C5]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR2C5
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	HIV
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein fragment corresponding to amino acids 1-231 of HIV P24 produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	25.58kDa



[View online »](#)