

## Product datasheet for **TA592304**

### **SARS-CoV-2 N Protein Rabbit Monoclonal Antibody [Clone ID: OTIR4E12]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR4E12
Applications:	ELISA, WB
Recommended Dilution:	WB 1:2000-5000, ELISA 1:5000-1:10000
Reactivity:	SARS-CoV-2
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Full length recombinant SARS-Cov-2 N protein produced in E.coli.
Formulation:	PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Database Link:	<a href="#">Entrez Gene 43740575 SARS-CoV-2</a>




[View online »](#)

Product images:



Western blot analysis of full length recombinant SARS-CoV-2 N protein (Cat# [TP790189], 0.02ug) by using anti-SARS-Cov-2 N protein antibody (TA592304).



Western blot analysis of full length recombinant SARS-CoV-2 N protein (Cat# [TP750180], 0.02ug) by using anti-SARS-Cov-2 N protein antibody (TA592304).