

## Product datasheet for **TA592330**

### SARS-CoV-2 Spike Protein Rabbit Monoclonal Antibody [Clone ID: OTIR2H1]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTIR2H1
Applications:	ELISA, WB
Recommended Dilution:	WB 1:2000-5000, ELISA 1:5000-10000
Reactivity:	SARS-CoV-2
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant fusion protein of SARS-CoV-2 Spike RBD Protein and Virion Surface Domain of Membrane Protein (M-protein) produced in HEK293
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Predicted Protein Size:	34kDa
Gene Name:	S Protein
Database Link:	<a href="#">Entrez Gene 43740568 SARS-CoV-2</a>




[View online »](#)

## Product images:



Western blot analysis of recombinant SARS-CoV-2 spike-RBD protein (Cat# [TP750182], 0.02ug) by using anti-SARS-CoV-2 spike protein antibody (Cat# TA592330).



Western blot analysis of recombinant fusion protein of SARS-CoV-2 spike RBD protein and virion surface domain of membrane protein (Cat# [TP701119], 0.05ug) by using anti-SARS-CoV-2 spike protein antibody (Cat# TA592330).