

Product datasheet for **TA592569AMS**

DDK Rabbit Monoclonal Antibody (Biotin conjugated) [Clone ID: OTIR5G2]

Product data:

| | |
|----------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | OTIR5G2 |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within DDK. The exact sequence is proprietary. |
| Formulation: | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. |
| Concentration: | 1mg/ml |
| Purification: | Purified from cell culture supernatant by affinity chromatography (protein A/G) |
| Conjugation: | Biotin |
| Storage: | Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws. |
| Stability: | Stable for 12 months from date of receipt |



[View online »](#)