

## Product datasheet for **TA592569JMS**

### DDK Rabbit Monoclonal Antibody (Magnetic Beads conjugated) [Clone ID: OTIR5G2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTIR5G2
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within DDK. The exact sequence is proprietary.
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	5mg/ml
Purification:	Purified from cell culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Magnetic Beads
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 12 months from date of receipt



[View online »](#)