

Product datasheet for TA592899

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Tau (MAPT) Rabbit Monoclonal Antibody [Clone ID: OTIR11C9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTIR11C9

Applications: ELISA

Recommended Dilution: ELISA 1:10000-1:100000

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Monoclonal

Immunogen: Synthetic acetylated peptide (the amino acid sequence is considered to be commercially

sensitive) within Human Tau (NP_005901). The exact sequence is proprietary.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

Predicted Protein Size: 45.7 kDa

Gene Name: microtubule associated protein tau

Database Link: NP 005901

Entrez Gene 4137 Human

P10636





Background: This gene encodes the microtubule-associated protein tau (MAPT) whose transcript

undergoes complex, regulated alternative splicing, giving rise to several mRNA species. MAPT transcripts are differentially expressed in the nervous system, depending on stage of neuronal maturation and neuron type. MAPT gene mutations have been associated with several neurodegenerative disorders such as Alzheimer's disease, Pick's disease,

frontotemporal dementia, cortico-basal degeneration and progressive supranuclear palsy.


[provided by RefSeq, Jul 2008]

Synonyms: DDPAC; FTDP-17; MAPTL; MSTD; MTBT1; MTBT2; PPND; PPP1R103; TAU; tau-40

Protein Families: Druggable Genome

Protein Pathways: Alzheimer's disease, MAPK signaling pathway

Product images:



ELISA with coating antigen Ac-tau274 synthetic acetylated peptide and 11C9 Ac-tau274 antibody (TA592899). Negative control: coating Full length tau protein