

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA592964

### PNPLA3 Rabbit Monoclonal Antibody [Clone ID: OTIR2E4]

### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR2E4
Applications:	IHC
Recommended Dilution:	IHC 1:2000
Reactivity:	Human
Host:	Rabbit
lsotype:	lgG
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human PNPLA3(NP_079501) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Predicted Protein Size:	52.8 kDa
Gene Name:	patatin like phospholipase domain containing 3
Database Link:	<u>NP_079501</u> <u>Entrez Gene 80339 Human</u> <u>Q9NST1</u>
Background:	The protein encoded by this gene is a triacylglycerol lipase that mediates triacylglycerol hydrolysis in adipocytes. The encoded protein, which appears to be membrane bound, may be involved in the balance of energy usage/storage in adipocytes. [provided by RefSeq, Jul 2008]
Synonyms:	ADPN; C22orf20; iPLA(2)epsilon




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

**Protein Pathways:** 

Glycerolipid metabolism, Glycerophospholipid metabolism, Limonene and pinene degradation, Metabolic pathways, Phenylalanine metabolism, Tyrosine metabolism

## **Product images:**



IHC staining of FFPE human kidney tissue within the normal limits using anti-PNPLA3 rabbit monoclonal antibody (TA592964) and Polink-2 HRP polymer detection kit ([D22-110]). Heatinduced epitope retrieval by EDTA solution buffer pH 8.0 ([B04C-100]) at 120°C for 3 min. The brown stain indicates positive stain, blue is the counter stain.(1:2000)

IHC staining of FFPE human liver tissue within the normal limits using anti-PNPLA3 rabbit monoclonal antibody (TA592964) and Polink-2 HRP polymer detection kit ([D22-110]). Heatinduced epitope retrieval by EDTA solution buffer pH 8.0 ([B04C-100]) at 120°C for 3 min. The brown stain indicates positive stain, blue is the counter stain.(1:2000)

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US