

Product datasheet for **TA600024**

IL8 (CXCL8) Mouse Monoclonal Capture Antibody [Clone ID: SNAP8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	SNAP8
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant protein of human IL-8 produced in E. coli.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	C-X-C motif chemokine ligand 8
Database Link:	NP_000575 Entrez Gene 3576 Human P10145
Synonyms:	GCP-1; GCP1; IL8; LECT; LUCT; LYNAP; MDNCF; MONAP; NAF; NAP-1; NAP1
Matched ELISA Pair:	TA700024
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Bladder cancer, Chemokine signaling pathway, Cytokine-cytokine receptor interaction, Epithelial cell signaling in Helicobacter pylori infection, NOD-like receptor signaling pathway, Pathways in cancer, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway




[View online »](#)

Product images:



IL-8 ELISA with SNAP8 Capture (TA600024) and 8CH Detection ([TA700024]) Antibodies. Substrate used: Recombinant Human IL-8 (E. Coli expressed Peptide)



IL-8 Luminex ELISA with SNAP8 Capture (TA600024) and 8CH Detection ([TA700024]) Antibodies. Substrate used: Recombinant Human IL-8 (E. Coli expressed Peptide)