

Product datasheet for **TA600033**

BIRC5 Mouse Monoclonal Capture Antibody [Clone ID: OTI3G3]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3G3
Applications:	ELISA
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human BIRC5 (NP_001159) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	baculoviral IAP repeat containing 5
Database Link:	NP_001159 Entrez Gene 332 Human O15392
Synonyms:	API4; EPR-1
Matched ELISA Pair:	TA700033
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Colorectal cancer, Pathways in cancer



[View online »](#)

Product images:



BIRC5 ELISA with 3G3 Capture (TA600033) and 2D9 Detection ([TA700033]) Antibodies. Substrate used: Recombinant Human BIRC5 ([TP305935])