

## Product datasheet for TA600037

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CDCP1 Mouse Monoclonal Capture Antibody [Clone ID: OTI2G2]

## **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2G2
Applications: ELISA

Recommended Dilution: ELISA 1:100-1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CDCP1 (NP\_073753) produced in HEK293T

cell

**Formulation:** Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose

**Concentration:** 0.5 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** CUB domain containing protein 1

Database Link: NP 073753

Entrez Gene 64866 Human

Q9H5V8

Synonyms: CD318; SIMA135; TRASK

Matched ELISA Pair: TA700037

**Protein Families:** ES Cell Differentiation/IPS, Transmembrane

