

## Product datasheet for **TA600050**

### DDK Mouse Monoclonal Capture Antibody [Clone ID: 1B12]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	1B12
Applications:	ELISA
Recommended Dilution:	ELISA 1:100-1:1000
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Matched ELISA Pair:	TA700098, TA700134, TA700135, TA700222, TA700223, TA700251, TA700252, TA700253, TA700254, TA700255, TA700256, TA700257, TA700258, TA700259, TA700260, TA700309, TA700310, TA700324, TA700325, TA700416, TA700486, TA700487, TA700488, TA700498, TA700514, TA700517, TA700518, TA700528, TA700535, TA700536



[View online »](#)