

## Product datasheet for TA600090

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### Factor XIIIa (F13A1) Mouse Monoclonal Capture Antibody [Clone ID: OTI3F1]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3F1

**Applications:** ELISA, LMNX **Recommended Dilution:** 1:100 - 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human F13A1 (NM\_000129) produced in HEK293T

cell

**Formulation:** Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose

**Concentration:** 0.5 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** coagulation factor XIII A chain

Database Link: NP 000120

Entrez Gene 60327 RatEntrez Gene 74145 MouseEntrez Gene 2162 Human

P00488

Synonyms: F13A

Matched ELISA Pair: TA700088


Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Complement and coagulation cascades






# **Product images:**



F13A1 ELISA with 3F1 Capture (TA600090) and 8E2 Detection ([TA700088]) Antibodies. Substrate used: Recombinant Human F13A1 ([TP306464])



F13A1 Luminex with 3F1 Capture (TA600090) and 8E2 Detection ([TA700088]) Antibodies. Substrate used: Recombinant Human F13A1 ([TP306464])