

Product datasheet for **TA600109**

LCMT1 Mouse Monoclonal Capture Antibody [Clone ID: OT11E7]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OT11E7
Applications:	LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human LCMT1 (NM_016309) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	leucine carboxyl methyltransferase 1
Database Link:	NP_057393 Entrez Gene 30949 Mouse Entrez Gene 361643 Rat Entrez Gene 51451 Human Q9UIC8
Synonyms:	CGI-68; LCMT; PPMT1
Matched ELISA Pair:	TA700109, TA700110
Protein Pathways:	Androgen and estrogen metabolism, Histidine metabolism, Selenoamino acid metabolism, Tyrosine metabolism




[View online »](#)

Product images:



LCMT1 Luminex with 1E7 Capture (TA600109) and 1G5 Detection ([TA700109]) Antibodies. Substrate used: Recombinant Human LCMT1 ([TP300018])



LCMT1 Luminex with 1E7 Capture (TA600109) and 4A5 Detection ([TA700110]) Antibodies. Substrate used: Recombinant Human LCMT1 ([TP300018])