

Product datasheet for **TA600124**


SOX17 Mouse Monoclonal Capture Antibody [Clone ID: OTI2G8]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI2G8
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 177-414 of human SOX17 (NP_071899) produced in E.coli.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SRY-box transcription factor 17
Database Link:	NP_071899 Entrez Gene 20671 Mouse Entrez Gene 312936 Rat Entrez Gene 64321 Human Q9H6I2
Synonyms:	VUR3
Matched ELISA Pair:	TA700123, TA700125
Protein Families:	Transcription Factors
Protein Pathways:	Wnt signaling pathway




[View online »](#)

Product images:


SOX17 ELISA with 2G8 Capture (TA600124) and 2F9 Detection ([TA700123]) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])



SOX17 Luminex with 2G8 Capture (TA600124) and 2F9 Detection ([TA700123]) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])



SOX17 Luminex with 2G8 Capture (TA600124) and 3H5 Detection ([TA700125]) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])