

## Product datasheet for **TA600147**

### ARH (LDLRAP1) Mouse Monoclonal Capture Antibody [Clone ID: OTI4B10]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI4B10
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human LDLRAP1 (NM_015627) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	low density lipoprotein receptor adaptor protein 1
Database Link:	<a href="#">NP_056442</a> <a href="#">Entrez Gene 100017</a> <a href="#">MouseEntrez Gene 500564</a> <a href="#">RatEntrez Gene 26119</a> <a href="#">Human Q5SW96</a>
Synonyms:	ARH; ARH1; ARH2; FHCB1; FHCB2
Matched ELISA Pair:	TA700147
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis




[View online »](#)

Product images:



LDLRAP1 Elisa with 4B10 Capture (TA600147) and 7G2 Detection ([TA700147]) Antibodies. Substrate used: Recombinant Human LDLRAP1 ([TP306643])



LDLRAP1 Luminex Elisa with 4B10 Capture (TA600147) and 7G2 Detection ([TA700147]) Antibodies. Substrate used: Recombinant Human LDLRAP1 ([TP306643])