

Product datasheet for **TA600183**

STAT4 Mouse Monoclonal Capture Antibody [Clone ID: OTI2A2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2A2
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human STAT4 (NM_003151) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	signal transducer and activator of transcription 4
Database Link:	NP_003142 Entrez Gene 20849 Mouse Entrez Gene 367264 Rat Entrez Gene 6775 Human Q14765
Synonyms:	SLEB11
Matched ELISA Pair:	TA700183, TA700192, TA700193, TA700194, TA700195
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Jak-STAT signaling pathway




[View online »](#)


Product images:




STAT4 ELISA with 2A2 Capture (TA600183) and 1B1 Detection ([TA700192]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])




STAT4 ELISA with 2A2 Capture (TA600183) and 1G9 Detection ([TA700183]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])




STAT4 ELISA with 2A2 Capture (TA600183) and 2F3 Detection ([TA700195]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])




STAT4 ELISA with 2A2 Capture (TA600183) and 2G5 Detection ([TA700193]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])




STAT4 ELISA with 2A2 Capture (TA600183) and 4H6 Detection ([TA700194]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])




STAT4 Luminex Elisa with 2A2 Capture (TA600183) and 1G9 Detection ([TA700183]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])




STAT4 Luminex Elisa with 2A2 Capture (TA600183) and 1B1 Detection ([TA700192]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])



STAT4 Luminex Elisa with 2A2 Capture (TA600183) and 2G5 Detection ([TA700193]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])



STAT4 Luminex Elisa with 2A2 Capture (TA600183) and 4H6 Detection ([TA700194]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])



STAT4 Luminex Elisa with 2A2 Capture (TA600183) and 2F3 Detection ([TA700195]) Antibodies. Substrate used: Recombinant Human STAT4 ([TP306892])