

## Product datasheet for **TA600228**

### **TYRO3 Mouse Monoclonal Capture Antibody [Clone ID: OTI2G2]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI2G2
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TYRO3 (NM_006293) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	TYRO3 protein tyrosine kinase
Database Link:	<a href="#">NP_006284</a> <a href="#">Entrez Gene 22174 Mouse</a> <a href="#">Entrez Gene 25232 Rat</a> <a href="#">Entrez Gene 7301 Human</a> <a href="#">Q06418</a>
Synonyms:	BYK; Dtk; Etk-2; Rek; RSE; Sky; Tif
Matched ELISA Pair:	TA700225, TA700227, TA700228, TA700229
Protein Families:	Druggable Genome, Protein Kinase




[View online »](#)


Product images:




TYRO3 ELISA with 2G2 Capture (TA600228) and 5B4 Detection ([TA700226]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])




TYRO3 ELISA with 2G2 Capture (TA600228) and 9B12 Detection ([TA700225]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])




TYRO3 Luminex Elisa with 2G2 Capture (TA600228) and 5B4 Detection ([TA700226]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])




TYRO3 Luminex Elisa with 2G2 Capture (TA600228) and 5B7 Detection ([TA700227]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])



TYRO3 Luminex Elisa with 2G2 Capture (TA600228) and 4C8 Detection ([TA700228]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])



TYRO3 Luminex Elisa with 2G2 Capture (TA600228) and 2D11 Detection ([TA700229]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])



TYRO3 Luminex Elisa with 2G2 Capture (TA600228) and 9B12 Detection ([TA700225]) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])