

Product datasheet for **TA600264**

Cystatin S (CST4) Mouse Monoclonal Capture Antibody [Clone ID: OTI1B1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1B1
Applications:	LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CST4 (NM_001899) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	cystatin S
Database Link:	NP_001890 Entrez Gene 1472 Human P01036
Synonyms:	MGC71923
Matched ELISA Pair:	TA700262



[View online »](#)

Product images:



CST4 Luminex Elisa with 1B1 Capture (TA600264) and 1A3 Detection ([TA700262]) Antibodies. Substrate used: Recombinant Human CST4 ([LY400707])