

# **Product datasheet for TA600268**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### TLE 1 (TLE1) Mouse Monoclonal Capture Antibody [Clone ID: OTI3H3]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3H3

Applications: ELISA, LMNX

Recommended Dilution: 1:100 - 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TLE1 (NM\_005077) produced in HEK293T

cell

**Formulation:** Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose

**Concentration:** 0.5 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** TLE family member 1, transcriptional corepressor

Database Link: NP 005068

Entrez Gene 21885 MouseEntrez Gene 362533 RatEntrez Gene 7088 Human

Q04724

**Synonyms:** ESG; ESG1; GRG1


**Matched ELISA Pair:** TA700268, TA700269, TA700270, TA700271

**Protein Families:** Transcription Factors







## **Product images:**



TLE1 Elisa with 3H3 Capture (TA600268) and 1F5 Detection ([TA700268]) Antibodies. Substrate used: Recombinant Human TLE1 ([TP304037])





TLE1 Elisa with 3H3 Capture (TA600268) and 2D8 Detection ([TA700269]) Antibodies. Substrate used: Recombinant Human TLE1 ([TP304037])




TLE1 Elisa with 3H3 Capture (TA600268) and 1D6 Detection ([TA700270]) Antibodies. Substrate used: Recombinant Human TLE1 ([TP304037])












TLE1 Luminex Elisa with 3H3 Capture (TA600268) and 1D6 Detection ([TA700270]) Antibodies. Substrate used: Recombinant Human TLE1 ([TP304037])




TLE1 Luminex Elisa with 3H3 Capture (TA600268) and 1F5 Detection ([TA700268]) Antibodies. Substrate used: Recombinant Human TLE1 ([TP304037])









TLE1 Luminex Elisa with 3H3 Capture (TA600268) and 3C9 Detection ([TA700271]) Antibodies. Substrate used: Recombinant Human TLE1 ([TP304037])